

How to Use Anatomage Share

Table of Contents:

What is Anatomage Share?	1
How to Create an Anatomage Share Account	1
How to Navigate the Home Screen	2
How to View, Like, Comment, Bookmark, and Report Files	4
How to Download Public Content from Anatomage Share	5
How to Upload Content to Anatomage Share	5
How to Edit your Uploaded Files	7
How to Delete an Uploaded File	7
How to Share a File with Another Anatomage Share User	8
How to Use Groups	9
How to Adjust your Account Details	11
How to Interact with User Pages	13
How to Download Presets from Anatomage Share onto the Table	14
How to Open a Downloaded File on the Table	16
How to Download Files into Anatomage Tablet or VR	17
How to Download Files into Anatomage Tablet or VR	17


What is Anatomage Share?

Anatomage Share is a website platform where users can access Anatomage Created materials (like example activities, tutorial videos, How-to Documents, and more), upload and store their Anatomage content, and share material with other Anatomage users around the world.

How to Create an Anatomage Share Account

Step 1: Go to <https://anatomageshare.com>


Step 2: Tap **Create Account**.

The image shows the Anatomage SHARE login and account creation page. At the top is the logo "Anatomage SHARE" with a double-headed arrow icon. Below the logo are two input fields: "E-mail" and "Password". The "Password" field has a small eye icon to its right. Below these fields is a "Log In" button. Under the "Log In" button are two links: "Remember Me" (with a checked checkbox) and "Forgot Password". Below these links is a paragraph of text about cookies: "We use cookies to ensure that we give you the best experience on our website. All cookies used are for site functionality and are deemed essential. By continuing to use our website, you agree with our use of cookies detailed in our [Privacy Policy](#)." Below this text is a horizontal line with the word "OR" in the center. At the bottom of the form, there is a link "New to Anatomage Share" followed by a link "Create Account". The "Create Account" link is circled in red.


Step 3: This will take you to an account creation page where you can input your name, email, institution information, and list what Anatomage products you use. Once you have filled in your information, agree to Anatomage's terms and tap **Create Account** to finish setting up your Anatomage Share account.

How to Navigate the Home Screen

1. Tap your name to see your profile and account information.
2. Tap the bell icon to see your notifications.
3. Tap Anatomage News to see recent Anatomage social media.
4. Tap Upload to upload content to Anatomage Share.
5. Tap Home to return to the Home Screen from another page.
6. Tap My Uploads to see all content you have uploaded (public and private)
7. Tap Shared with me to see all content other users have shared with you.
8. Tap Groups to see Groups you are a part of, or join and create new groups.
9. Tap Bookmarks to see all content you have bookmarked.
10. Tap Follow List to see the list of AnatomageShare users that you follow.
11. Tap Contact Us to submit questions and feedback to Anatomage.
12. Tap Help to see FAQs.




13. Select a filter tab to filter to only resources that fit that category.
14. Type in the Search Bar to search through available resources by keyword or creator name. Tap the Advanced Search button to add filters like tags and file location.
15. Tap Anatomage Top 10 arrow to view all the resources in the current Top 10.
16. Tap on any resource to open in.
17. Tap on an Uploader name to see their personal page.



How to View, Like, Comment, Bookmark, and Report Files

Step 1: Open the file you'd like to view details for by clicking on it.

1. Scroll through preview images of the file using the arrows.
2. View the file name, description, uploader and date.
3. Tap the Like button to like a file
4. Tap the 3 dots to report a file to Anatomage.
5. Tap the bookmark icon to save a resource to your Bookmarks.
6. Type a comment and tap post to comment on the file.




How to Download Public Content from Anatomage Share


Step 1: Locate the resource you'd like to download, and tap on it to open it.

Step 2: If you would like to download all files in the folder, or if there is only one document to download, tap "Download All" to begin the download automatically.


Note: You may need to allow pop ups on anatomageshare.com in order to download multiple files at once!




Step 3: If you would only like to download select files from a grouped folder, hover over the file you would like until a download icon pops up. Select this icon to begin the download.







How to Upload Content to Anatomage Share

Step 1: Tap  **Upload** to begin the upload process.


Step 2: Hover over **Supported File Formats**  to see a list of all file formats supported by Anatomage Share.

Step 3:

- Drag and drop files to the  **Drag & Drop Files** box to drop them directly into Anatomage Share.
- Tap  **Browse Files** to upload a single file from your computer.
- Tap  **Browse Folder** to upload a full folder of files from your computer.
- Tap  **Video Link** to embed a web video directly into AnatomageShare.

Step 4: Give your file or folder a Title and a detailed Description.

Step 5: Choose which Tag(s) apply to your files.



Step 6: Decide if you want your folder to be private (just for you) or public (**available to all of Anatomage Share**) by using the Visibility check box. Publicly shared files will not count toward your free private storage limit.



Note: Files are listed as public by default.

Step 7: Optional - upload a Thumbnail as your file/folder's preview image.



Step 8: Edit the Author Name if needed.




Step 9: Once all required fields are filled in, tap Upload



Step 10: Go to My Uploads to see all uploaded files.

How to Edit your Uploaded Files

Step 1: Tap My Uploads to see all your uploaded files. Click on a particular file to open the Content Detail View for that file.

Step 2: To change the visibility preference of a file (public or private) from the My Uploads menu - tap the current visibility status and choose the one you would like.



Step 3: To change other details of a file from the My Uploads Menu - hover over the file and tap the  icon, then choose  Edit Content Details .

To change other details of a file from the Content Detail View - tap the edit button in the top right




Step 4: From here you can adjust all file settings, including Title, Description, Tags, Visibility, Thumbnail, and Author.

Step 5: Tap Save to save any changes to the file.

How to Delete an Uploaded File

Step 1: Go to the my Uploads menu.


Step 2: Tap the dot  next to the file you would like to remove.

Step 3: Tap  Remove in the top right corner.

How to Share a File with Another Anatomage Share User

From the My Uploads page:

Step 1: Tap My Uploads to see your uploaded file.

Step 2: Hover over the file you would like to send and tap the Airplane Icon .

Step 3: Enter the email address(es) you would like to share the file to.

Step 4: Tap Add to add the email address.

Step 5: Tap Send to send the file. The recipient will receive an email and a notification in their Anatomage Share account if they have one. The file will appear in their “Shared With Me” menu.

From the Content Detail View:

Step 1: Tap My Uploads to see your uploaded file.

Step 2: Click on the File you would like to share. You can share public and private files. The file will go directly to the “Shared with me” folder of the recipient.


Step 3: Tap the Airplane icon to share the file.




Step 4: Follow Steps 3-5 above.

How to Use Groups

How to Make a Group

Step 1: Open the Group menu in the left toolbar.  Group

Step 2: Under “Create Group”, enter a Name and Description for your Group.




Group Name*

Sports Medicine


Group is created as private and users the owner's storage.

Group Description*

Lab and lecture material for High School Sports Medicine courses. Focus on Mu:

 Create Group

Step 3: Tap Create Group. This will automatically generate a group code. **Share this code** with other users so they can join a group.




Group: Sports Medicine

Join Code

4RTNPHEF


Copy

ⓘ If you reset the code, the current code will not be valid to join the group.


 [Code Reset](#)

Close

How to Join a Group

Step 1: Open the Group menu in the left toolbar.  Group

Step 2: Under “Join a Group with a Code”, enter the group code that corresponds to the group you would like to join (this must be provided to you by another group member).



Step 3: Tap Join Group to join.

How to Upload to a Group

Upload Directly to the Group



Step 1: In Group, open the page for the Group you would like to upload to.

Step 2: Tap  Group Upload

Step 3: Upload content as normal (see How to Upload Content to Anatomage Share).

Add a Public Document to the Group


Step 1: Go to the detail page for the resource you would like to add to the group.

Step 2: Tap the three dots  and choose Add to Group. 

Step 3: Select the Group you would like to link the resource to from the list.

Add a Private Document to the Group

Step 1: Go to My Uploads to see your available resources.

Step 2: Hover over the resource you would like to add, and tap the three dots , and choose Add to Group


OR: Go to the Content Detail View for the resource, and follow the same steps as adding a public document to the group.

How to Edit or Delete a Group

Step 1: In Group, open the page for the Group you would like to edit or delete (must be your group).


Step 2: Tap the Settings icon. 

Step 2: Tap Edit Group Information to edit the Group. Tap Delete Group to delete the Group.




How to Leave a Group

Step 1: In Group, open the page for the Group you would like to leave.

Step 2: Hover over the Joined button. 

Step 3: Tap Leave.



How to Adjust your Account Details

How to Edit your Information

Step 1. Tap your name in the top right corner.

Step 2: Select My Account.

Step 3: Tap the  icon to make edits.


Step 4: Tap the  icon to save edits.

How to Adjust your Home Tags

Step 1. Tap your name in the top right corner.

Step 2: Select My Account

Step 3: Tap the Home Tags drop down and select any tags you would like to prioritize on your Home Page.




How to Adjust your Notification Preferences

Step 1. Tap your name in the top right corner.

Step 2: Select My Account

Step 3: Toggle your notification preferences.



How to Change your Password

Step 1. Tap your name in the top right corner.

Step 2: Select My Account.

Step 3: Enter your current password and select Change Password to reset.

How to Delete your Anatomage Share Account

Step 1. Tap your name in the top right corner.

Step 2: Select My Account

Step 3: Select Delete Account - Note: All uploaded public and private files, shared files, and comments will also be removed. All data is non-recoverable.

How to Interact with User Pages

How to View your Page


Step 1. Tap your name in the top right corner.

Step 2: Select My Page

How to View Someone Else's Page

Step 1: Tap on the Name of the person whose Page you want to view.

Names can be found on the preview of content and in the content detail view.



Skeletal System Quiz - Bones of the Skull
15 question quiz of the major b...

Jordan Sipla
Anatomage - Others

Jordan Sipla
Anatomage - Others

Skeletal System Quiz - Bones of the Skull

0 Downloaded | Updated: 12/29/2022

0 Likes

How to Follow a Page

Step 1: Go the Page of the person you want to follow

Step 2: Tap Follow

+ Follow

How to Download Presets from Anatomage Share onto the Table

Presets You Uploaded



Step 1: Upload the Preset to AnatomageShare.com following the instructions in this Document.

Step 2: Connect your Anatomage Table to the Internet via Ethernet (available on all Tables) or WiFi (Available on Tables shipped after April 2021).

Step 3: In the main application Menu, choose Anatomage Share



Step 4: Log-in with the same ID and Password you use for the Anatomage Share website.



Step 5: Go to Files

Step 6: Go to the "My Uploads" tab.

Step 7: Check the box next to a VPF file to download that file to the Table



Presets Others Uploaded

Step 1: On your AnatomageShare.com account, mark the VPF you want to access on the Table in one of three ways:


1. Bookmark it with the Bookmarks Tab.
2. Have the file shared within a Group you are in.
3. Have the file creator share the file with you directly.

Step 2: Connect your Anatomage Table to the Internet via Ethernet (available on all Tables) or WiFi (Available on Tables shipped after April 2021).

Step 3: In the main application Menu, choose Anatomage Share



Step 4: Log-in with the same ID and Password you use for the Anatomage Share website.

A dark gray dialog box titled "Anatomage Share" with a close button (X) in the top right corner. It contains two input fields: "Email:" and "Password:", each with a small "X" icon to its right. Below the fields is a paragraph of text: "Connect your Anatomage Table online and sign in to your Anatomage Share account. Anatomage Share allows easy storage and sharing of preset and DICOM files between the Anatomage Table, Anatomage VR, and other Anatomage products. To create an account, visit www.anatomageshare.com." Below this is another paragraph: "All content downloaded from Anatomage Share is for informational and educational purposes only and does not constitute medical advice." At the bottom is an "Okay" button.

Step 5: Go to Files

Step 6: Depending on how the file was shared to you, go to the "Bookmarks", "Shared", or "Groups" tab.

Step 7: Check the box next to a VPF file to download that file to the Table



How to Open a Downloaded File on the Table

Through the Anatomage Share Menu


Step 1: Follow the previous instructions to download the VPF file onto the Table.

Step 2: In the main application Menu, choose Anatomage Share



Step 3: Log-in with the same ID and Password you use for the Anatomage Share website.

Step 4: Go to Files



Step 5: Go to the tab in which the VPF file is located



Step 6: Tap on the title of the VPF file you would like to open. This will automatically open the scan the VPF was created in, and pre-load the presets into the preset menu.

Through the Preset Menu

Step 1: Open the Cadaver or Scan in which you would like to open a downloaded preset.


Step 2: Open the Preset settings menu



Step 3: Open the Folders icon



Step 4: Choose the Anatomage Share folder



Step 5: Tap the preset you would like to open to open it on the Table.

Note: All presets downloaded onto the Table will appear here. If you attempt to open a preset designed in a different cadaver than the one you have opened, the Table will warn you, but proceed with opening the file.

How to Download Files into Anatomage Tablet or VR

(For Anatomage Tablet / VR customers only)

Anatomage VR must be connected to the internet in order to access Anatomage Share.




Step 1: In Anatomage VR main menu, tap the Account icon.

Step 2: Log-in with the same ID and Password you use for the Anatomage Share website.

Step 3: Tap Files to access uploaded files.



Step 4: Tap the download icon to download files onto Anatomage VR.



The download icon will turn gray when a file has been fully downloaded to AnatomageVR.

How to Download Files into Anatomage Tablet or VR

(For Anatomage Tablet / VR customers only)

Files must be downloaded to Anatomage VR before they can be opened. See previous step for instructions.

For .VPF Preset Files:

Step 1: Open the cadaver or Case Library scan that you would like to open the preset in.





Step 2: Tap the Preset Menu.


Step 3: Tap the Import icon  to see available VPFs.

Step 4: Locate your preset in the Import Preset Menu under Anatomage Share Files.


For Practice Quiz .VPF Files:

Step 1: Open the cadaver (male or female) that you would like to Practice in.  

Step 2: Select the Practice icon 

Step 3: Tap the import icon .

Step 4: Locate your preset in the Load Quiz Menu under Anatomage Share Files.



For DICOM Files:

Step 1: Locate the downloaded DICOM series in the files menu.  → Files

Step 2: Double tap the DICOM file you downloaded. It should open automatically.